

SomnoBrux®

The most comfortable device for the treatment of Bruxism



- **QUALITY** – superior quality design, manufacturing and materials
- **COMFORT** – offers maximum comfort and durability yet discreet to wear and does not limit tongue movement
- **CUSTOM-FIT** – perfectly tailored to suit the individual on either top or bottom arch
- **SMH BFLEX SOFT LINER** – provides a comfortable fit with no odour or discolouration
- **ONE YEAR WARRANTY** – on materials to give you peace of mind

For advice or to order, please contact:



SomnoBrux®

Like all SomnoMed products, the SomnoBrux benefits from superior quality design, manufacturing and materials. Custom made by our experienced SomnoMed® dental technicians to fit to the anatomy of the patient using SMH BFlex soft liner.

SMH BFlex Lining

The lining of the SomnoBrux® is soft and comfortable thanks to the use of our unique resin SMH BFlex which also ensures excellent retention. The SomnoBrux® device covers all occlusal surfaces of the dental arch and ensures conformity of contact. The nature of the lining prevents odour and discolouration.

Our patented SMH BFlex is a flexible material with permanent elasticity which forms a chemical and mechanical bond with the rigid and flexible materials of the device. This ensures an exceptional level of fit precision, retention, comfort and durability.

Fitting

We recommend providing an impression of each arch and a bite registration in order to provide the necessary occlusal thickness for the acrylic resin. This makes it possible to achieve a perfectly tailor-made device and ensures better comfort for your patient. It's you who decides which arch should wear the device.



SomnoMed AG is ISO13485-certified

Somnomed products are registered in:

EU, Switzerland, United States, Canada, Australia, South Korea,
Japan, Hong Kong

Somnomed products are patented medical devices.

www.somnomed.com