

SomnoMed Signature Line

Clinically proven world class oral devices



Avant™

CAD-CAM manufacturing with BFlex soft liner and innovative titration system that has interchangeable straps providing an 8 mm range and adjustments at 1 mm increments.




Herbst Advance Elite™

Medicare/PDAC approved, BFlex liner with 18% more Flexural Strength. 8 mm range of calibration with SomnoMed's proprietary visual calibration indicator.




Fusion®

Rapid advancement with interchangeable wings providing an 8.5 mm range and fine screw adjustments at .1 mm increments.



Flex®

Features a soft inner liner made with SMH BFlex material, a proprietary soft polymer that is molecularly bonded to the acrylic.



Classic®

SomnoMed's original customized winged oral device, the mandible is held in place using medical grade stainless steel ball clasps.

To order, call or email SomnoMed Sales
888.447.6673, option 4, ussales@somnomed.com

 SomnoMed®


SomnoMed Standard Line

Treatment focused. Technology driven.



Herbst Advance®

Medicare/PDAC approved, 8 mm range of calibration with SomnoMed's proprietary visual calibration indicator.



Air/Air+

Air available with ball clasps or standard Poly+ liner for retention. Air+ available with standard Poly+ liner for retention. Affordably priced, ideal for lower reimbursement scenarios.



Morning Repositioner

Hands-free bite deprogrammer designed to help return the patient's mandible back to its pre-treatment centric position. Can be ordered with any device and is fabricated off original impressions or scans.

To order, call or email SomnoMed Sales
888.447.6673, option 4, ussales@somnomed.com