

SomnoDent[®] HERBST ADVANCE ELITE[™]

BFlex Soft Inner Liner

- Superior comfort.



18% Greater Flexural Strength from Precured PMMA

3 Year “No Questions Asked” Replacement Plan

- If it breaks, we replace it.

PDAC (Medicare E0486) Approved

- The only milled, PDAC approved oral appliance with a soft inner liner.



What is Flexural Strength?

Flexural strength is the resistance to bending and twisting forces.


18% Greater Flexural Strength from Precured PMMA



Flexible, Strong & Durable

- PMMA flexes during movement.
- Lab tests show that SomnoMed's pre-cured PMMA is 18% stronger when subjected to flexing forces.


The SomnoDent HERBST ADVANCE ELITE is made of Precured PMMA.



Polymer/Monomer (Standard Acrylic)



Precured PMMA



Flexural or Bending Strength

Reference: SOM-REP-0649A, Table 7.1



3 Year “No Questions Asked” Replacement Plan

SomnoMed[®] makes it easy. Get a replacement SomnoDent HERBST ADVANCE ELITE appliance at no charge for 3 years from date of manufacture if your patient's device breaks.



www.somnomed.com | 903506 Rev A